


Stirring Rods, paddle type


Simax Borosilicate Glass Catalogue Number **3100**



DIAMETER MM	LENGTH MM	PACK CONTENT
6	125	10
6	150	10
6	200	10
6	300	10

Plain Rod, polished ends


Simax Borosilicate Glass Catalogue Number **7765**



DIAMETER MM	LENGTH MM	PACK CONTENT
6	100	10
6	200	10
6	250	10
6	300	10
6	400	10

Test tubes with rim


Borosilicate Glass Catalogue Number **7801**



d MM	L MM	PACK CONTENT
10	75	100
10	100	100
12	75	100
12	100	100
16	100	100
16	125	100
16	150	100
18	150	100
24	150	50

Test tubes with rim, medium wall


Simax Borosilicate Glass Catalogue Number **7830**



d MM	L MM	WALL MM	PACK CONTENT
10	75	1.0	100
10	100	1.0	100
12	75	1.0	100
12	100	1.0	100
16	100	1.0	100
16	125	1.2	100
16	150	1.2	100
18	150	1.2	100
24	150	1.2	50

Culture tubes without rim


Borosilicate Glass Catalogue Number **3800**



d MM	L MM	PACK CONTENT
10	75	1000
12	75	1000
12	100	1000
16	125	1000
16	150	1000

Test tubes without rim, medium wall


Simax Borosilicate Glass Catalogue Number **3840**



d MM	L MM	WALL MM	PACK CONTENT
10	75	1.0	100
12	75	1.0	100
12	100	1.0	100
16	125	1.0	100
16	150	1.2	100
24	150	1.2	50

Test tubes with screw cap


Borosilicate Glass Catalogue Number **3801**



d MM	L MM	WALL MM	PACK CONTENT
16	100	1.8	100
16	120	1.8	100
16	160	1.8	100

Test tubes with screw cap


Simax Borosilicate Glass Catalogue Number **3830**



GL	d MM	L MM	WALL MM	PACK CONTENT
18	16	100	1.8	100
18	16	120	1.8	100
18	16	160	1.8	100

Tubes centrifuge, graduated


Simax Borosilicate Glass Catalogue Number **3822**



CAPACITY ML	d MM	L MM	PACK CONTENT
10	17	120	50
15	18	140	50
25	24	150	50

Test tubes with conical bottom, with rim

Simax Borosilicate Glass Catalogue Number **3821**



d MM	L MM	WALL MM	PACK CONTENT
14	100	1.5	50
14	120	1.5	50
15	80	1.2	50
15	105	1.2	50
16	100	1.8	50
16	120	1.8	50
17	120	1.2	50